

# nice Neotech Medical Systems Pvt. Ltd.

## nice 5050 Respiration Monitor



### OPERATING/INSTALLATION MANUAL

This operating manual provides all the information necessary for the user to safely set up and operate this equipment. It is the responsibility of the user to follow the instructions and recommendations provided

Part Number - 73-00-040

Revision - 01

Date - 01/02/2017



85. Krishna Industrial Estate, Vannagaram, Mettukuppam Chennai-95.

Ph: 91-44-24764608, Telefax: 91-44-24762594

E-mail: [info@niceneotech.com](mailto:info@niceneotech.com) / [marketing@niceneotech.com](mailto:marketing@niceneotech.com) / [service@niceneotech.com](mailto:service@niceneotech.com)

Web: [www.niceneotech.com](http://www.niceneotech.com)

## **Contents**

**User Responsibility/Operator Profile**

**Declaration for Language**

**Model Description**

**Definitions**

**Definition of Warning Indication**

**Section A: Warnings**

**Section B: Cautions**

**Section C: Symbols**

**Section 1: Description**

1.1 Intended Use

1.2 General

1.3 Control unit

**Section 2: Installation**

2.1 Setup

2.2 Installation of Respiration Sensor

2.3 Charging of Monitor

2.4 Audio and Visual Indication System Test

**Section 3: Operation**

3.1 Control Panel Operation

3.2 Attaching the Respiration Sensor

3.3 Monitor Setup

3.4 Silencing the 'clicking sound'

3.5 Changing the Audio and Visual Indication delay

3.6 Checking Audio and Visual Indication output

3.7 Monitoring the infant

3.8 Switching off the Monitor

3.9 If the Audio and Visual Indication sounds

**Section 4: Cleaning and Maintenance**

4.1 Cleaning

4.2 Lifetime of the product

4.3 Do's & Don'ts

**Section 5: Specification**

**Section 6: Warranty**

**Section 7: Trouble Shooting**

7.1 System Failure

7.2 Maintenance interval

7.3 Disposing of the equipment

**Section 8: Spare Parts List**

## **User Responsibility/Operator profile**

This Product will perform in conformity with the description thereof contained in this operating manual and accompanying labels and/or inserts, when assembled, operated, maintained and repaired in accordance with the instructions provided. This Product must be checked periodically. A defective Product should not be used. Parts that are broken, missing, plainly worn, distorted or contaminated should be replaced immediately. Should such repair or replacement become necessary, nice Neotech recommends that a telephone or written request for service advice be made to the nearest nice Neotech Regional Service Center.

This Product or any of its parts should not be repaired other than in accordance with written instructions provided by nice Neotech and by nice Neotech trained personnel. The Product must not be altered without the prior written approval of nice Neotech's Quality Assurance Department. The user of this Product shall have the sole responsibility for any malfunction which results from improper use, faulty maintenance, improper repair, damage, or alteration by anyone other than nice Neotech.

**Warning:** Before using the nice Neotech Monitor, read this entire manual. Attempting to use this device without a thorough understanding of its operation may result in patient or user injury. This device should only be operated by personnel trained in its operation and under the direction of qualified medical personnel familiar with the benefits and risks of this type of device.

## **Declaration for Languages**

User Manual and label will be provided in the appropriate language to ensure that the user understands. Language validation will be done for the language of the user manual, Label, Corresponding documents, when nice Neotech Medical Systems Private Limited supplies to EU countries.

## **Model Descriptions:**

nice 5050 Respiration Monitor.

Definitions:

**Select:** choose the Audio and Visual Indication setting

**Self Test:** Booting the equipment to check the internal peripheral.

**APNEA:** is a term for suspension of external breathing

**APNEA Period:** The delay time in seconds for each breath.

**Respiration:** The Process of the inspiration & exhalation of gas from the body.

**Definition of Warning indication:**

Three levels of warning indication are used throughout this manual and on the unit. They are defined as follows,

A **DANGER** notice indicates an immediately hazardous situation which, if not avoided, will result in death or serious injury, serious damage to property such as total loss of use of equipment, and a fire.

A **WARNING** notice indicates an indirectly (Potentially) hazardous situation which if not avoided, will result in death or Serious injury, serious damage to property such as total loss of use of equipment, and a fire.

A **CAUTION** notice indicates a hazardous situation which, if not avoided can result in minor or moderate injury, partial damage to property and loss of data stored in computers.

## Section A: Warnings

**Warning:** The Monitor should only be used under the direction of a clinician. Failure to observe this warning may cause infant injury or death.

**Warning:** The Monitor must be operated within the specified parameters, e.g. temperature, see Technical Specification section 5. Failure to observe this warning could affect the functions of the Monitor and may cause infant injury or death.

**Warning:** The performance of the Monitor is dependent on the selection of the correct Monitoring rate for the infant and the environment of use. Failure to observe this warning could affect the functions of the Monitor and may cause infant injury or death.

**Warning:** Excessive movement may adversely affect the Monitor's performance. Care must be taken when used in an ambulatory situation and the Monitor placed in a secure position free from excessive movement. Failure to observe this warning could affect the functions of the Monitor and may cause infant injury or death.

**Warning:** Before using the nice Neotech Respiration Monitor, read this entire manual. Attempting to use this device without a thorough understanding of its operation may result in patient or user injury.

**Warning:** Do not use the monitor in the presence of flammable anesthetics; a possible explosion hazard exists under these conditions.

**Warning:** In a dust free area, keep hands clean and then install the equipment.

**Warning:** Do not perform the Checkout Procedures (Mechanical and Control Unit) with a patient.

**Warning:** Complete the "Checkout Procedures" section of this manual before putting the unit into operation. If the Monitor fails any portion of the checkout procedures it must be removed from use and repaired.

**Warning:** There is a possibility of electromagnetic interference or other interference causes from other external equipment, the monitor may get in operation. Use EMC compliance Equipment to avoid the interference

**Warning:** The Monitor may cause radio interference, in which case adequate measure may be required to prevent interference

**Warning:** Mobile communication devices e.g. mobile phones, should not be used at a distance of less than one metre to the Monitor, as this could affect the safe operation of the Monitor. Even if these devices cause no obvious interference, possible interference can only be ruled out if this advice is strictly observed. Malfunction of the Monitor can cause incorrect or loss of monitoring resulting in infant injury or death

**Warning:** The Monitor is not suitable for use on infants with pacers. Failure to observe this warning may affect the functions of the Monitor and may cause infant injury or death.

**Warning:** The Monitor is not suitable for use with adults or children over the age of 18 months. Failure to observe this warning may cause injury or death.

**Warning:** The Monitor must not be covered by bedding or other objects so that audible and visual Audio and Visual Indication are impaired. Failure to observe this warning could result in infant injury or death.

**Warning:** To prevent the Monitor being accidentally switched off, it must not be placed within reach of the infant. Failure to observe this warning could result in infant injury or death.

**Warning:** The Monitor must not be pushed into confined spaces, possibly resulting in pressure being exerted on the keypad buttons causing it to be switched off or the settings changed. Failure to observe this warning could result in infant injury or death.

**Warning:** If the main audible and visual Audio and Visual Indication are activated, immediately attend to the infant. This is of prime importance before any other actions are taken. In the event of any Audio and Visual Indication, the infant must first be made safe and comfortable. Failure to observe this warning may cause infant injury or death.

**Warning:** Use only Rechargeable batteries. alkaline batteries are unsuitable and must not be used as they

could affect safe operation of the Monitor. Failure to observe this warning could affect the functions of the Monitor and may cause infant or user injury.

**Warning:** When the Low Battery indicator is lit, it is essential to ensure re-charge. Failure to observe this warning could affect the functions of the Monitor and may cause infant or user injury.

**Warning:** Following a significant liquid spill onto the Monitor, it should be wiped dry and inspected by an authorized technician before being returned to clinical use. Failure to observe this warning could affect the functions of the Monitor and may cause infant or user injury.

**Warning:** The Respiration Sensor should only be used on one infant and for a maximum of 7 days. Using the same sensor on more than one infant increases the risk of cross infection. Failure to observe this warning may cause infant injury or death

## Section B: Cautions

**Caution:** Isolation from the supply mains is Detachable charger provided

**Caution:** Use cleaning solution sparingly on a cloth when cleaning the Monitor. Do not saturate the unit - excessive solution causes damage to internal components.

**Caution:** Do not autoclave or gas sterilize the Monitor & sensor. Do not immerse the device & sensor in liquid cleaner.

**Caution:** Only competent individuals trained in the repair of this equipment should attempt to service it as detailed in the service manual. The Service Manual provides detailed information solely for use by individuals having proper knowledge, tools and test equipment, and for service representatives trained by nice Neotech.

**Caution:** The equipment may affect while using the defebulator.

**Caution: Use of nonstandard components:** Consult the manufacturer for repair and replacement of components. Use of incorrect component can adversely affect Safety, performance and/or damage the equipment performance.

**Caution:** Clean surface using wet cloth dipped in mild soap water and squeeze dry excess water before use.

**Caution:** Do not allow water to spill into equipment.

**Caution:** Avoid using any solvent, spirits, alcohol to clean plastic parts.















**Caution:** Do not use water or any other liquid to clean Electronics and Electrical parts.














**Caution:** Do not autoclave, steam sterilise, ETO sterilise or subject the Monitor to temperatures in excess of 45° C (113° F). Excessive temperatures may cause damage to the Monitor.

**Caution:** Do not use unapproved cleaning agents.

**Caution:** Do not lift the Monitor by the sensor tubing.

## Section C: Symbols

Mark	Indication
	Caution
	Caution, electric shock
	Caution, ignition (Indicates possible ignition under certain condition)
	Warning
	Warning, Do not touch
	Warning, Don't use fire
	Protective Earth (ground)
	Alternating current
	Off (Power: disconnection from main)
	On (Power: connection to the main)
	"Off" (only for a part of Equipment)
	"On" (only for a part of Equipment)
	Audio Indicator
	Audio Indication Silence

	Serial Number
	Authorised Representative in European Community
	Date of Manufacture
	Manufacturer
	Catalogue Number
	Refer Instruction for use
	Drip proof protection to IPX1
	Type B Equipment
	Mute Key
	Select Key
	ON/OFF Key
	Charging Indication
	Low Battery Indication

## **Section 1: Description**

1.1 Intended Use

1.2 General

1.3 Control unit

### **1.1 Intended Use**

Respiration Monitor is intended to continuously monitoring the respiration of infants up to the age of 18 months and will give audio visual Audio and Visual Indication if the infant stops breathing for more than a pre-set time.

### **1.2 General**

When an infant is being monitored it is essential that the user remains within visual and audible range of the Monitor so that the Audio and Visual Indication can be seen and heard.

### **1.3 Control Unit**

The breathing signal is picked up by a soft, foam-filled Sensor (Respiration Sensor) that is attached to the infant's abdomen. During breathing, expansion of the abdominal wall causes a small amount of air to pass through the Sensor and tube to the Monitor. This triggers an audible click and simultaneous visual indication.

## Section 2: Installation

### 2.1 Setup

### 2.2 Installation of Respiration Sensor

### 2.3 Charging of Monitor

### 2.4 Audio and Visual Indication System Test

## 2.1 Setup

After removal from the shipping containers, inspect the monitor and all accessory items for any signs of damage which may have occurred during shipment. File a damage claim with the shipping carrier if damage has occurred. Also confirm the presence of all accessory items or factory installed options as listed on the packing slip.

## 2.2 Installation of Respiration Sensor



Picture 1

- Plug one end of the respiration sensor in to the device sensor port.

## 2.3 Charging of Monitor



Picture 3

- For the initial charging keep the device for charging continuously 12 hours up to the green LED flashes

## 2.4 Audio and Visual Indication system test


- ❖ The operation of the Audio and Visual Indication system must always be checked as follows before the Monitor is used:
- ❖ Prior to connecting a Respiration Sensor, press the ON/OFF key to on the device.
- ❖ Select 10 sec by using select key. The monitor confirms the selection of apnea period by a short beep and the ON 10 sec light briefly illuminates.
- ❖ Do not apply any pressure on the sensor pad
- ❖ The audio visual Audio and Visual Indication (flashing red Audio and Visual Indication indicator) should trigger after 10 seconds.
- ❖ The Monitor is turned off by pressing the HOLD FOR OFF button. (Keep the button pressed until you hear a short beep).

❖ If the Audio and Visual Indication system test is unsuccessful, see section 7 trouble shooting guide.

### Section 3: Operation

- 3.1 Control Panel Operation
- 3.2 Attaching the Respiration Sensor
- 3.3 Monitor Setup
- 3.4 Silencing the 'clicking sound'
- 3.5 Changing the Audio and Visual Indication delay
- 3.6 Checking Audio and Visual Indication output
- 3.7 Monitoring the infant
- 3.8 Switching off the Monitor
- 3.9 If the Audio and Visual Indication sounds

#### 3.1 Control Panel Operation (See Picture 4)



Picture 4


#### Front Panel Button and Indicators

1	Audio and Visual Indication Indication	7	Charging Indication
2	Mute On	8	Battery Low Indication
3	10 Second Indication	9	ON/OFF Key
4	Select Key	10	Sensor Port

5	15 Second Indication	11	Buzzer
6	20 Second Indication		

### 3.2 Attaching the Respiration Sensor

- Clean and dry the infant's abdomen.
- The most suitable position for attaching the Respiration Sensor is just below the navel, slightly to the left or right. Fix the Respiration Sensor pad and tube using hypo allergenic tape 50 mm wide and 100 mm long. Position the Sensor tube down between the infant's legs to the Monitor, so that it is away from the neck, or so that the infant cannot play with it (see illustration). The Respiration Sensor can be worn under clothing (even nappies). In this case, however, ensure that the Sensor tube is passed through the clothing.
- It is important to ensure that the Respiration Sensor connector is pushed firmly into the Monitor's connector socket using a slight twisting action. Do not apply excessive pressure.



Picture 1

- **NOTE:** When attaching the Respiration Sensor pad ensure that the tape is long enough and does not interfere with the Respiration Sensor. The tape should not be stretched too tightly over the Sensor pad. Only the Neotech recommended Respiration Sensor should be used.

### 3.3 Monitor set up

- Switching on the Monitor
- With the Respiration Sensor attached to the Monitor, switch the Monitor on by pressing the ON/OFF key, and select the required APNEA Period (10, 15 & 20 sec Audio and Visual Indication delay) by using select key.
- The choice of which Audio and Visual Indication delay should always be determined by a clinician
- The ON 10 sec. Audio and Visual Indication delay, will Audio and Visual Indication if no breath is detected after 10 seconds.
- The ON 15 sec. Audio and Visual Indication delay, will Audio and Visual Indication if no breath is detected after 15 seconds
- The ON 20 sec. Audio and Visual Indication delay, will Audio and Visual Indication if no breath is detected after 20 seconds.
- For every breath detected, the Monitor will make a clicking sound accompanied by a green ON 10 sec or a Yellow ON 15 sec or an orange ON 20 sec flashing light.
- When the Monitor is first switched on there may be a small delay before it starts to detect the infant's breaths.

- If the selected Audio and Visual Indication time is breached an audiovisual Audio and Visual Indication is triggered, this can be silenced by pressing mute key.
- Should the infant start to breathe again before mute key is pressed then the audible Audio and Visual Indication will stop, but the red Audio and Visual Indication light will continue to flash for another 60 seconds.

### 3.4 Silencing the 'clicking sound'

The clicking which sounds each time a breath is detected, can be silenced by simultaneously pressing the Mute key and the select key.

To restore the click repeat the above procedure. However, breathing continues to be indicated by the flashing indicator lamp

### 3.5 Changing the Audio and Visual Indication delay

If the Monitor is set to an Audio and Visual Indication delay interval of 10 seconds, 15 seconds or 20 second Audio and Visual Indication interval is required. This can be selected by pressing the select key.

### 3.6 Checking Audio and Visual Indication output

To check the Audio and Visual Indication output, remove the Respiration Sensor connector from the Monitor; the audible Audio and Visual Indication should sound after 10 seconds (15 seconds or 20 seconds), and the Audio and Visual Indication indicator flashes. Replace the Respiration Sensor connector, the audible Audio and Visual Indication stops and the visual Audio and Visual Indication continues to flash for about 60 seconds. The monitor immediately reverts to the monitoring function.

However, in order to save battery power it is recommended that the visual Audio and Visual Indication be cancelled as soon as possible by pressing the mute key

### 3.7 Monitoring the infant

While the infant is being monitored by the respiration monitor it is essential that the user remains within audible and visual range of the Monitor. If using the Monitor in an infant's cot position it away from the infant i.e. at the bottom of the cot and hang the Monitor from the cot bars. The infant can also be monitored at any other time, for example: in a buggy, car seat etc.

### 3.8 Switching off the Monitor

To switch the Monitor off press the ON/OFF key (for about 3 seconds) until a short beep is heard.

### 3.9 If the Audio and Visual Indication sounds

**Warning:** If the main audible and visual Audio and Visual Indication are activated, immediately attend to the infant. This is of prime importance before any other actions are taken. In the event of any Audio and Visual Indication, the infant must first be made safe and comfortable. Failure to observe this warning may cause infant injury or death.

- Go to the infant but do not panic.
- If the infant has stopped breathing carry out the resuscitation techniques that you have learned and immediately seek medical assistance.
- If the red Audio and Visual Indication indicator is flashing but the audible Audio and Visual Indication is silent, this shows that a previous Audio and Visual Indication condition was present. Press mute key to cancel the visual Audio and Visual Indication indication and observe the infant for a while. If an audible Audio and Visual Indication sounds again, carry out the Audio and Visual Indication system test see Section 2.4, Audio and Visual Indication system test.
- If the audible Audio and Visual Indication is repeated, press mute key again and observe the infant.

- If the infant is safe and breathing normally but the Monitor is not indicating the infant's breath, there is possibly a technical fault, see section 7 "Trouble Shooting".

## Section 4: Cleaning & Maintenance

### 4.1 Cleaning the Monitor

### 4.2 Lifetime of the product

### 4.3 Do's & Don'ts

**Warning:** Periodically check the insulation and the connection of the charger cable it may cause fire because of poor insulation and short circuits due to aging.

**Warning:** Switch off the equipment and disconnect the Power cord from the mains before take in to cleaning

**Warning:** Don't pour the water for cleaning it may enter into the electronics circuits it cause short circuit and get shock.

**Caution:** Use the cleaning solution sparingly on a cloth when cleaning the Monitor. Do not saturate the unit - excessive solution causes damage to internal components.

**Caution: Use of nonstandard components:** Consult the manufacturer for repair and replacement of components. Use of incorrect component can adversely affect Safety, performance and/or damage the equipment performance

The following lists recommended cleaning solutions that may be used safely.

Generic Formulation	Max Concentration level
Hydrogen peroxide	6%
Sodium hypochlorite compound	0.5 % Aqueous Solution
Quaternary ammonium	0.2 %
Iodophors	50 mg/liter
Glutaraldehyde	2%

**Caution:** Use of cleaning/disinfecting solutions containing chemicals that are not listed above (i.e. alcohol, acetone, etc.), or chemicals in greater concentrations than those listed above, may damage the material being cleaned.

**Caution:** Do not autoclave or gas sterilize Monitor and Respiration sensor.

The instruction below describes the methods to clean and disinfect the instrument and its accessories.

### 4.1 Cleaning the Monitor

Cleaning- The external surfaces of the Monitor body may be cleaned by wiping them with a cloth moistened with a mild detergent solution.

### 4.2 Lifetime of product

Since the product is classified under programmable medical electrical system and incase of unavailability of microcontroller the life time of the product can be considered as minimum five years.

### 4.3 Do's and Don'ts

#### DO

- Read all of the directions before using for the first time.
- Visually inspect the sensor for damage before each use.
- Test batteries regularly and charge the device when battery indicator shows low battery
- Keep the unit, sensor and connections are securely tight
- Clean the case with isopropyl alcohol or mild detergent only

## **DON'T**

- Use this Monitor if you suspect any malfunction.
- Use the device in the presence of flammable gases.
- Use alkaline batteries.
- Autoclave or freeze the sensor or device.
- Immerse the unit or sensor in any liquid.
- Expose the unit to devices that produce high levels of radio, short wave, microwave, x-ray, or high frequency interference.
- Use cleaning agents or liquids in the cable receptacles or around the battery compartment
- Place the unit itself in a water vapor-saturated environment.
- Expose the unit to a condensing water environment such as a mist tent.

## Section 5: Specifications

Apnea interval	: 10, 15 & 20 seconds
Sensor type	: Passive Pneumatic
Audio Indication	: Breath & Apnea Audio and Visual Indication
Visual Indication	: Breath & Apnea Audio and Visual Indication
Battery Type	: 7.2VDC Li-Ion Rechargeable
Battery Indication	: Low Battery & Charging
Apnea Audio and Visual Indication Indication sec (Orange)	: 10 sec (Green), 15 sec (Amber) & 20
Power Supply	: Battery / Main Supply
Battery Charger	: 10.5VDC/1.5A
Battery Low	: Visual indication will activate at low battery
Battery backup	: 24 hours continuous
Charging Time	: 7 hours continuous
Dimensions	: 9cm (L) X 6cm (D) X 17.5 cm (H)
Weight	: 550 grams

## Environment

Operating Temperature	: 15°C to 35°C
Operating Relative Humidity	: 15 – 90% RH, non condensing
Storage Ambient Temperature	: -10°C to 60°C
Storage Relative Humidity	: 50 - 90% RH, non-condensing
Altitude	: Sea level to 1.9 miles 3 Km

## IEC 60601-1 Specifications


Type of Protection against electric Shock	: Class 1
Degree of Protection against electric Shock	: Type BF
Mode of Operation	: Continuous
Protection against hazardous of explosion	: Not protected
Protection against increase of liquid	: IPX1
Product Classification	: Class IIb.

**Handling and Storage:**



nice Neotech Medical Systems Pvt Ltd  
No 85, Krishna Industrial Estate  
Vannagaram, Mettukuppam  
Chennai - 600095, INDIA

**HANDLE WITH CARE**  
**LIFE SAVING MEDICAL EQUIPMENT**



## Section 6: Warranty

### 6.1 Conditions

1. The warranty is confined to the first purchaser of the product only and is not transferrable.
2. Repairs under warranty period shall be carried out by the company authorized personnel only
3. In the event of repairs of any part/s of the unit, this warranty will thereafter continue and remain in force only for the unexpired period of the warranty. The time taken for repair and in transit whether under the warranty or otherwise shall not be excluded from the warranty period.
4. In case of any damage to the product/misuse detected by the Authorized service personnel the warranty conditions are not applicable and repairs will be done subject to availability of parts and on a chargeable basis only
5. Wear and Tear, and defects caused by manipulation or unsuitable treatment are not included under the warranty.
6. We warranty this unit for 12 months from the date of Installation. Warranty includes the repair and replacement of faulty components.
7. Defects caused by improper use, and defects due to causes beyond control like lightning, abnormal voltage, acts of god, and also defects caused by rats, cockroaches or any other insects will not be covered under warranty.
8. Warranty is not applicable if the equipment is not purchased from Neotech/authorized Neotech Dealer
9. Warranty is not applicable if the warranty card is not filled and sent back to Neotech.
10. Equipment has an expected shelf life of 5 years and service life of 6 years.

### Customer Details cum Warranty Card

Date: \_\_\_\_\_

Hospital Name & Address: \_\_\_\_\_

\_\_\_\_\_

Contact Person & Telephone/Fax No \_\_\_\_\_

Email \_\_\_\_\_

Department: NICU / PICU / OT / Gynecology / Causality / Others \_\_\_\_\_

Equipment Name: \_\_\_\_\_

Model No: \_\_\_\_\_ Sl. No. \_\_\_\_\_

Date of Purchase: \_\_\_\_\_ Date of Installation \_\_\_\_\_

Name of Authorized Dealer: \_\_\_\_\_

Customer Signature & Date  
(I accept the terms & conditions of Warranty)

Dealer Signature with seal

-----

Kindly fill the above and sent the same

From \_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

To:  
The Service In-charge  
nice Neotech Medical Systems Pvt. Ltd.  
No, 85. Krishna Industrial Estate,  
Vannagaram, Mettukuppam,  
Chennai-600095. Tamil Nadu,INDIA  
Ph: 91-44-24762594, 24764608  
Tele Fax: 91-44-24766920 – Ext: 203  
Email: [service@niceneotech.com](mailto:service@niceneotech.com), [info@niceneotech.com](mailto:info@niceneotech.com)  
Web: [www.niceneotech.com](http://www.niceneotech.com)  
Toll Free No. 1800-425-2594 (India only)

## Section 7: Trouble Shooting

7.1 System Failure

7.2 Maintenance interval

7.3 Disposing of the equipment.

### 7.1 System Failure

Sl. No	Problem	Cause	Remedy
1	Monitor does not work.	No batteries fitted.  Batteries fitted incorrectly.  Batteries discharged.	Fit new batteries.  Check that the batteries have been fitted as indicated in the battery compartment.  Recharge / Fit new batteries.
2	No respiration signals.	Respiration Sensor not connected to the Monitor.  Sensor pad has become loose on the skin.	Make sure that the Sensor connector is inserted in the Monitor's connector socket, using a twisting motion.  Remove the Sensor; clean and dry the infant's skin. Attach the Sensor pad with new hypo allergenic tape. See Section 3.2, attaching the Respiration Sensor.
3	Monitor has become defective.	Internal fault.  Respiration Sensor is defective.	Test by removing the Sensor tube from the Monitor (with a slight turning movement) and by pressing gently with a finger over the inlet socket. This should produce signals on the Monitor. If it does not, the Monitor is faulty and must not be used. Contact nice Neotech  Attach a new Respiration Sensor to the Monitor, see Section 2.2, Attaching the Respiration Sensor. If the Respiration Sensor still does not work, the Sensor may be faulty. Contact nice Neotech.
4	No display	Batteries expired.  Bad battery connection.  Charger defects	Check/recharge the battery.  Check battery connections.  Replace the charger
5	Keys inoperable/cannot turn unit off	Electronics component failure	Remove the rear panel & disconnect the battery

Note: In the event that none of these procedures produce desired results, remove the batteries and return the unit to nice Neotech for repair

### 7.2 Maintenance Intervals:

- Always disinfect and clean the unit and accessories before any maintenance – even when returning the unit to the supplier for repair.
- Always disconnect power supply before any maintenance.

- Use only nice Neotech's original parts for maintenance.

**Warning:** Periodically check the insulation and the connection of the charger cable it may cause fire because of poor insulation and short circuits due to aging.

**Warning:** Don't misaligned the EMI Shielding and the beads it may cause the EMI interference to the equipment

### **7.3 Disposing of the Unit**

- At the end of its Service life Dispose of the equipment in accordance with National waste Disposal Regulations or ask a suitable Disposal contractor to dispose of the unit. The local Environmental agency can supply further details.

## Section 8: Spare Parts List

Sl. no	Part No	Part Name	Qty	Unit
1	50-05-053	PCB Assembly	1	No.
2	88-00-107	Baby Respiration Sensor	1	No.
3	91-00-137	Battery	1	No.
4	91-00-138	Charger	1	No.

### Service contact:

#### nice Neotech Medical Systems Pvt. Ltd.

85. Krishna Industrial Estate, Vannagaram,  
Mettukuppam Chennai-600095. Tamil Nadu, INDIA.

Ph: 91-44-2476 4608 Telefax: 91-44-2476 2594

E-mail:[service@niceneotech.com](mailto:service@niceneotech.com) / [info@niceneotech.com](mailto:info@niceneotech.com)

Web:[www.niceneotech.com](http://www.niceneotech.com) /

[www.cemark.info/mdd/niceneotech.html](http://www.cemark.info/mdd/niceneotech.html)

### EU Authorised Representative

CMC Medical Devices

C/Horacio Lengo N° 18

P 29006, Malaga – Spain

Mr. Manuel Mateos

Manager – Regulatory Affairs

Tel: +34951214054

Fax: +34952330100

[mmateos@cmcmmedicaldevices.com](mailto:mmateos@cmcmmedicaldevices.com)

[www.cmcmmedicaldevices.com](http://www.cmcmmedicaldevices.com)

**For Complaints/Comments/Suggestions or Feedbacks**

Date: \_\_\_\_\_

Hospital Name & Address: \_\_\_\_\_

\_\_\_\_\_

Contact Person & Telephone/Fax No. \_\_\_\_\_

Email \_\_\_\_\_

Department: NICU / PICU / OT / Gynecology / Causality / Others \_\_\_\_\_

Equipment Name: \_\_\_\_\_

Model No: \_\_\_\_\_ Sl. No. \_\_\_\_\_

Date of Purchase: \_\_\_\_\_ Date of Installation \_\_\_\_\_

Complaints \_\_\_\_\_

\_\_\_\_\_

Comments: \_\_\_\_\_

\_\_\_\_\_

Suggestions: \_\_\_\_\_

\_\_\_\_\_

Feed back: \_\_\_\_\_

\_\_\_\_\_

-----  
Kindly fill the above and sent the same

From:

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

To:

The Marketing In-charge

nice Neotech Medical Systems Pvt. Ltd.

No, 85. Krishna Industrial Estate,

Vannagaram, Mettukuppam,

Chennai-600095. Tamil Nadu,

INDIA.

Ph: 91-44-24762594, 24764608

Tele Fax: 91-44-24766920 – Ext: 203

Email: [marketing@niceneotech.com](mailto:marketing@niceneotech.com)

Toll Free No. 1800-425-2594 (India only)